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Introduction

The **tCAP(N3)** reagent is designed to modify antibodies by conjugating the **Ac-Lys(N₃)-Gly-Gly** peptide in a buffer at pH 7–9. This reagent specifically introduces two molecules of Ac-Lys(N₃)-Gly-Gly into the **Lys248 residue** of the antibody's **Fc region** (one molecule per heavy chain). [1] The introduced azide groups enable the preparation of antibody-drug conjugates (ADCs) or fluorescently labeled antibodies via downstream click chemistry. Structurally, the tCAP(N3) reagent consists of three functional components (**Fig1.(i)**). The preparation and evaluation of ADCs using tCAP chemistry have also been reported [2].

1. **An azide tripeptide (Ac-Lys(N₃)-Gly-Gly)** that serves as the payload modified onto the antibody.
2. An active acyl donor site.
3. **An Fc-specific affinity peptide moiety.**


Figure 1. Modification Mechanism of the tCAP(N3) Reagent

When the tCAP(N3) reagent is added to an antibody solution in a neutral buffer, the affinity peptide directs and positions the reagent near the Fc region. Subsequently, the side-chain amino group of the antibody's lysine (Lys248) residue nucleophilically attacks the active site. This reaction transfers the azide tripeptide onto the antibody, yielding the azide-functionalized antibody (IgG-N₃). Because the reaction occurs at the active acyl donor site located between the affinity peptide and the payload within

the reagent structure, only the azide tripeptide payload is successfully introduced into the antibody (**Fig1. (ii)**). Although the affinity peptide is released after the transfer reaction, their affinity for the antibody is maintained, so they remain attached to the antibody via non-covalently bound. Therefore, to completely eliminate the affinity peptide from the system, it must be dissociated using a dissociation buffer at approximately pH 3.3, followed by removal techniques such as dialysis, ultrafiltration, or gel filtration column. **Even if some affinity peptides remain, they are attached only via reversible noncovalent bonds; therefore, they are not expected to significantly affect the antibody's binding to FcRn or its pharmacokinetics.**

Protocols

[Reagent and Sample Preparation]

- IgG Concentration: Prepare the IgG solution in PBS at a concentration of 1–10 mg/mL (approx. 6.7–67 μ M)
- 100 μ M tCAP(N3) in H₂O stock solution
(tCAP(N3) Stock Solution (100 μ M): Dissolve in H₂O according to the vial's net weight (tCAP(N3) MW = 4798.2) to obtain a 100 μ M concentration. Example: Add 229 μ L of H₂O for a 0.11 mg vial)

Caution: Do not dissolve the tCAP(N3) reagent in a neutral buffer. The reagent will decompose and become inactive.

Storage: Store the 100 μ M tCAP(N3) stock solution in H₂O in a -20°C freezer.

- 200 mM carbonate-bicarbonate buffer (pH 8.9)

[Materials for Affinity Peptide Removal]

- 10% of AcOH aq.
- 150 mM NaCl/0.1%AcOH aq.
- dialysis tube or a gel filtration column (such as a Sephadex® G-50 column).

Protocol 1: Standard Method (DAR: Max 2, but typically to afford 1.5–1.8)

1. **Mixing:** Add tCAP(N3) (100 μ M in H₂O) to the IgG solution at a **1:4–5 molar ratio** (antibody : reagent), and mix vigorously.
2. **Affinity Binding:** Incubate the mixture for 10 min to allow affinity binding.
Note: The addition of NaCl can be effective for antibodies that are prone to aggregation.
3. **pH Adjustment:** Add 200 mM carbonate-bicarbonate buffer (pH 8.9), relative to the initial IgG solution volume to adjust the pH.
Note: Please check the pH. A low pH will decrease the DAR.
4. **Conjugation:** Incubate the solution for 2 h at room temperature.
5. **Purification:** Depending on your requirements, please choose one of the following three

methods:

5-1. For Acid-Sensitive IgG:

If the target IgG is unstable under acidic conditions, perform the dialysis against a neutral pH buffer.

Note: Under neutral conditions, the affinity peptide cannot be removed, **leaving 2.0 equivalents** in the sample.

5-2. Standard Purification (Partial peptide removal)

Quench the reaction by adding 10 vol% of 10% AcOH aq. relative to the initial IgG solution volume (this adjusts the mixture to approximately pH 3.3). Dialyze the sample against a buffer containing 150 mM NaCl/0.1%AcOH aq. to remove unreacted peptides. Subsequently, perform a second dialysis against PBS to neutralize the pH.

Note: This standard method results in incomplete affinity peptide removal, leaving approximately **0.4 equivalents** of the peptide in the sample.

5-3. Advanced Purification (complete peptide removal)

Quench the reaction by adding 10 vol% of aqueous 10% AcOH aq. relative to the initial IgG solution volume (this adjusts the mixture to approximately pH 3.3). Purify the sample using a gel filtration column (such as a Sephadex® G-50 column, equilibrated and eluted with 150 mM NaCl/0.1%AcOH aq.), and collect the flow-through fraction. Subsequently, dialyze the collected fraction against a neutral buffer.

Note: This standard method achieves **complete affinity peptide removal**.

6. Yield: This process yields IgG-N₃ with a DAR of 1.5–1.8.

Technical Notes:

*To obtain IgG-N₃ with a higher drug-to-antibody ratio (**DAR ≈ 2.0**), the reaction must be repeated.*

Please contact us at (sales@peptide.co.jp) to obtain the detailed protocols.

Protocol 2: Downstream Click Reaction

4. **Mixing:** Mix the obtained IgG-N₃ with the DBCO-reagent at a 1:3 molar ratio.
5. **Conjugation:** Incubate the reaction mixture for 3~4 h at room temperature. (Extended reaction time may be required depending on the molecular size and structure of the DBCO compound.)
6. **Purification:** Remove the remaining unreacted DBCO reagents via dialysis or ultrafiltration.

Examples: Modification to Trastuzumab (Tmab)

[samples]

- Trastuzumab (Tmab) 0.5 mg Concentration : 8.7 mg/mL(= 58 μM) in PBS

- tCAP(N3) : 100 μ M tCAP(N3) in H₂O (add 229 μ L H₂O to 0.11 mg vial)
 - 200 mM carbonate-bicarbonate buffer (pH 8.9)
 - ↓ 8.7 mg/mL (= 58 μ M) of Trastuzumab (Tmab) in PBS = 58 μ L (3.36 nmol)
 - ↓ Add 100 μ M tCAP(N3) in H₂O = 134 μ L (13.4 nmol, 4 eq)
 - ↓ Stir 10 min at room temperature
 - ↓ Add 200 mM carbonate-bicarbonate buffer (pH 8.9) = 63 μ L
 - ↓ Stir 2 hours at room temperature
- [Final IgG concentration : 13 μ M (2.0 mg/mL), Final carbonate-bicarbonate buffer concentration 50 mM]
- ↓ Dialysis against neutral buffer
 - ↓ Check the DAR by LC-MS

IgG-N₃


Figure 2. LC-TOF-MS analysis results. Left) LC chromatogram detected at 280 nm. Right) MS spectrum of the IgG conjugate. While the retention time of T-Mab remained unchanged after the modification, the MS spectrum showed a distinct shift. This result indicated that approximately 1.8 equivalents of Ac-Lys(N₃)-Gly-Gly were successfully conjugated to the antibody, confirming a DAR of 1.8.

Product Information

Supported Antibody Species

Applicable Species: Human IgG1, IgG2, IgG4; Mouse IgG1, IgG2a, IgG2b, IgG3; Rabbit IgG, etc.

Note: The antibody affinity peptide of the tCAP(N3) reagent is derived from a Protein A fragment peptide. Therefore, any antibody capable of binding to Protein A is considered

compatible with this modification method.

Reaction pH and DAR

While the modification reaction can proceed at pH 7.0 or 8.0, performing the reaction under these conditions will result in a lower drug-to-antibody ratio (DAR). For optimal DAR results, following the recommended protocol pH is highly advised.

Notice of Usage Restrictions & Licensing

This product is sold under a patent license agreement with Kagoshima University. Utilization for radioactive pharmaceuticals is prohibited under this agreement.

Reference

- [1] H. Kawakami, A. Rafique, S. Tsuda, A. Nagashima, Y. Nohara, N. Dozono, K. Sakamoto, S. Masuda, M. Kiyoshi, H. Shibata, A. Ishii-Watabe, T. Yoshiya, and Y. Ito, *Chem.Eur.J.*, **32**, e03144 (2026).
- [2] K. Tsuji, Y. Miura, T. Kuwata, R. Matsuzaki, T. Kobayakawa, K. Matsumoto, Y. Ito, T. Yoshiya, S. Matsushita, and H. Tamamura, *ChemMedChem*, **21**, e202500820 (2026).